

Ocular Ahmed DVX Surgical Gonio Lens



The Ocular Ahmed Direct View eXtreme wide field lens is the latest redirecting surgical gonioscope available and it is exclusively from Ocular Instruments. Excellent for MIGS procedure. Widest double mirror field of view with 360° rotation and outstanding image clarity. Steam Autoclavable.

Designed for upright direct viewing WITHOUT TILTING the microscope.

*Designed with Iqbal (Ike) K. Ahmed, MD, FRCSC
Toronto, ON Canada
U.S. Patent #10,716,470*

Creating Solutions to Save Sight!

For more information on this product or our complete product line contact our

Customer Service Department or visit our website.

Take advantage of our 30 day no risk trial!



Comparison views of the Ocular Ahmed DVX and Ocular Swan Jacob Autoclavable Surgical Gonio lenses

Ocular Ahmed DVX Surgical Gonio

With the Ocular Ahmed DVX Surgical Gonio lens, it is NOT necessary for the physician to tilt the microscope.




Image represents a model eye



Incorrect view...Lens is tilted



CORRECT VIEW

Ocular Swan Jacob Autoclavable Surgical Gonio

With the Ocular Swan Jacob Autoclavable Surgical Gonio lens, it is necessary for the physician to tilt the microscope.




Image represents a model eye