

SecureFlex® HF Surgical Gonio Lens



"The first **Hands Free** gonioscopy lens for anterior chamber angle surgical procedures".

Exclusively from Ocular Instruments, Inc.

Made of Medical Grade Silicone, this flexible optic design resists decoupling from patients' eyes during surgical procedures.

Unique self retaining feature eliminates hand-held manipulation required by other competitive lenses; *creating a bi-manual surgical experience.*


Product Code: **OSIG**

Mag: **1.2x**

FOV: **90°**



Photo represents image of a model eye.



- Ideal for anterior chamber angle surgeries
- Visualization for MIGS type procedures
- Hands Free self retaining design
- Resists "decoupling" of optical interface during surgery
- Wide access to clear cornea
- Allows for Bimanual surgery
- Material - Medical Grade Silicone
- Flexible lens flange for stability
- Single use device
- Sold in a carton of 10

Available as a convenient single use device eliminating the risk of cross patient contamination.

Designed in collaboration with Ivantis Inc., Irvine, CA

U.S. Patent: 10,413,178



Creating solutions to save sight!



For more information on this product or our complete product line contact our Customer Service Department or visit our website.

